Android SDK API

Version: 0.8

Date: 2023/10/2



Category

Revision History	2
Introduction	
Specification	3
Connected to the Device	5
Perform Scan – Normal	9
Perform Scan – Quick Set	
Perform Scan – Manual	
Update Reference to the Device	
Activation Key Setting	
Device Information	17
Device Status	17
List Scan Configuration	
Add Scan Configuration to the Device	
Warm-Up the Device	19
Calculate AU from a New Reference Scan	20



Revision History

Revision	Author	Date	Description		
0.1	Iris	2020/5/26	1. Initial version		
0.2	Iris	2020/7/21	 Add lamp ramp up ADC and lamp ADC among repeated times. 		
0.3	Iris	2020/9/16	 Modify the flow chart of the Connected to the Device. Add the section of other receiver. 		
0.4	Iris	2021/4/22	 Add the flow chart of the all section. Modify the description of the all section. Remove the section of Lamp Ramp Up ADC and Lamp ADC among repeated times. Add to the section of Perform Scan – Normal, Perform Scan – Quick Set and Perform Scan – Manual. Remove the section of other receiver. Add to the section of Perform Scan – Manual. 		
0.5	Iris	2021/7/2	 Add the section of Warm-Up the device. Modify the flow chart of the Connected to the Device and add 1-(b) to teach how to do the warm-up. 		
0.6	Iris	2022/1/4	 Add the section of introduction and specification. Add note to the required section. Add corresponding API table in each section. 		
0.7	Iris	2023/3/31	 The specification items add extend plus machines. The Connected to the Device chapter removes the need to use Notify_IsEXTVersion() to notify the ISC SDK whether the device is an extended wavelength device. 		
0.8	Iris	2023/10/2	 Add Calculate AU from a New Reference Scan topic. 		



Introduction

This document is to let users understand how to use the SDK provided by ISC for development. It introduces the entire APP process, including connection, scanning..., and provides how to call API for development. Corresponding with source code will make it easier to get started. If you want to know more about BLE callback, please refer to ISC NIRScan Bluetooth Communications Data Sheet V1.3.pdf.

Specification

The TIVA versions and functions supported by this SDK are as follows:

TIVA	Function									
	Normal	Quick	Manual	Maintain	Activation	Device	Device	Read Device	Add Device	Lock
	Scan	Set Scan	Scan			Information	Status	Configuration	Configuration	Button
<v2.1.0.67< td=""><td>х</td><td>х</td><td>х</td><td>Х</td><td>x</td><td>x</td><td>х</td><td>х</td><td>х</td><td>х</td></v2.1.0.67<>	х	х	х	Х	x	x	х	х	х	х
v2.1.0.67~v2.3.2	V	х	х	Х	х	V	V	V	х	х
v2.4.2	V	V	V	V	V	V	V	V	V	Х
≧v2.4.3	V	V	V	V	V	V	V	V	V	V

Standard Wavelength (900nm~1700nm)

Extend Wavelength (1350nm~2150nm)

TIVA		Function								
	Normal	Quick	Manual	Maintain	Activation	Device	Device	Read Device	Add Device	Lock
	Scan	Set	Scan			Information	Status	Configuration	Configuration	Button
		Scan								
<v3.3.0< td=""><td>V</td><td>v</td><td>V</td><td>V</td><td>V</td><td>V</td><td>V</td><td>V</td><td>V</td><td>х</td></v3.3.0<>	V	v	V	V	V	V	V	V	V	х
≧v3.3.0	V	V	V	V	V	V	V	V	V	V



TIVA	Function									
	Normal	Quick	Manual	Maintain	Activation	Device	Device	Read Device	Add Device	Lock
	Scan	Set	Scan			Information	Status	Configuration	Configuration	Button
		Scan								
≧v5.0.1	V	V	V	V	V	V	V	V	V	V

Extend Plus Wavelength (1600nm~2400nm)



Connected to the Device

ScanViewActivity.java



- 1. When the device is connected, it will first check whether the user want to warm up the device.
 - (a) Call the ISCNIRScanSDK.SetCurrentTime() to synchronize the time. After synchronize the time, will download calibration coefficient and calibration matrix automatically. Register
 RefCoeffDataProgressReceiver to see the progress of download calibration coefficient. Register
 CalMatrixDataProgressReceiver to see the progress of download calibration matrix. Register
 RefDataReadyReceiver to get the calibration coefficient and matrix.

Note: In our SDK, should pass calibration coefficient and calibration matrix parameter to ISCNIRScanSDK library. Refer to the code under RefDataReadyReceiver :

refCoeff = intent.getByteArrayExtra(ISCNIRScanSDK.EXTRA_REF_COEF_DATA); refMatrix = intent.getByteArrayExtra(ISCNIRScanSDK.EXTRA_REF_MATRIX_DATA); ArrayList<ISCNIRScanSDK.ReferenceCalibration> refCal = new ArrayList<>(); refCal.add(new ISCNIRScanSDK.ReferenceCalibration(refCoeff, refMatrix)); ISCNIRScanSDK.ReferenceCalibration.writeRefCalFile(mContext, refCal);



API	Process	Receiver	Receiver
			Descriptio
			n
ISCNIRScanSDK.SetCurrentTime() (UUID: 0x4348410C-444C-5020-4E49- 52204E616E6F	Synchronize the time->Download calibration	RefCoeffDataProgressReceiver (UUID: 0x43484110-444C-5020-4E49- 52204E616E6F)	Download calibration coefficient
0x4348410F-444C-5020-4E49-52204E616E6F 0x43484111-444C-5020-4E49-52204E616E6F)	coefficient- >Download calibration matrix	CalMatrixDataProgressReceiver (UUID: 0x43484112-444C-5020-4E49- 52204E616E6F)	Download calibration matrix
		RefDataReadyReceiver	Get the calibration coefficient and matrix

(b) Call the ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.ON) to open the lamp to warm up the device then set Lamp_Info to LampInfo.WarmDevice. Register ReturnSetLampReceiver to do the next step.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK .LampState.ON) (UUID: 0x43484144-444C-5020-4E49- 52204E616E6F)	Open the lamp to warm up the device	ReturnSetLampRecei ver	Notify the completion of the lamp setting, according to lamp info to do next

2. Call the **ISCNIRScanSDK.GetActiveConfig()** to request to get active configuration from the device. Register **GetActiveScanConfReceiver** to get the index of active configuration.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetActiveConfig() (UUID: 0x43484118-444C-5020-4E49- 52204E616E6F)	Request to get active configuration	GetActiveScanConfReceiver (UUID: 0x43484118-444C-5020- 4E49-52204E616E6F)	Get the index of active configuration

 Call the ISCNIRScanSDK.GetScanConfig() to request to get the number of scan configuration and data from the device. Register ScanConfSizeReceiver to get the number of scan configuration. Register ScanConfReceiver to get the scan configuration.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetScanConfig() (UUID:	Get number of scan configuration ->Get	ScanConfSizeReceiver	Get the number of scan configuration
0x43484114-444C-5020-4E49-52204E616E6F)	data	ScanConfReceiver (UUID: 0x43484115-444C-5020-4E49- 52204E616E6F)	Get the scan configuration



4. Call the **ISCNIRScanSDK.GetSpectrumCoef()** to request to get the spectrum calibration coefficient. Register **SpectrumCalCoefficientsReadyReceiver** to get the spectrum calibration coefficient.

	API	Process	Receiver	Receiver Description
ſ	ISCNIRScanSDK.GetSpectrumCoef()	Request to get the	SpectrumCalCoefficientsReadyReceiver	Get the spectrum
	(UUID: 0x4348410D-444C-5020-	spectrum calibration	(UUID: 0x4348410E-444C-5020-4E49-	calibration coefficient
	4E49-52204E616E6F)	coefficient	52204E616E6F)	

Note: In our APP, should pass spectrum calibration coefficient to AddScanConfigViewActivity.java:

passSpectrumCalCoefficients = SpectrumCalCoefficients;

5. Call the ISCNIRScanSDK.GetDeviceInfo() to request to get the device information including the device of model name, serial number, HW version and TIVA version. Register DeviceInfoReceiver to get the device information. It is necessary to determine the device is standard, extend or extend plus wavelength according to TIVA version. If the device is a standard wavelength, the FW level is obtained with GetFWLevelStandard(Tivarev). If the device is a extension wavelength, the FW level is obtained with GetFWLevelEXT(Tivarev). If the device is a extension plus wavelength, the FW level is obtained with GetFWLevelEXT(Tivarev). Should call InitParameter() to initialize the maximum and minimum wavelengths.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetDeviceInfo() (UUID: 0x00002A29-0000-1000-8000- 00805F9B34FB 0x00002A24-0000-1000-8000-00805F9B34FB 0x00002A25-0000-1000-8000-00805F9B34FB 0x00002A27-0000-1000-8000-00805F9B34FB 0x00002A28-0000-1000-8000-00805F9B34FB	Request to get the device information	DeviceInfoReceiver	Get the device information
Notify_IsEXTVersion()	Notifying the ISC SDK that device is a standard wavelength or an extended wavelength	X	X
GetFWLevelStandard(Tivarev)	Get standard wavelength FW level	X	x
GetFWLevelEXT(Tivarev)	Get extend wavelength FW level	X	x
InitParameter()	Initialize the maximum and minimum wavelengths	X	x

If the TIVA version \geq 2.1.0.67 for the device, continue with the following steps.

6. Call the **ISCNIRScanSDK.GetMFGNumber()** to request to get the Manufacturing Serial Number of the device. Register **ReturnMFGNumReceiver** to get the Manufacturing Serial Number.



API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetMFGNumber() (UUID: 0x4348410B-444C-5020- 4E49-52204E616E6F)	Request to get the Manufacturing Serial Number	ReturnMFGNumReceiver (UUID: 0x4348410B-444C-5020- 4E49-52204E616E6F)	Get the Manufacturing Serial Number

7. Call the **ISCNIRScanSDK.GetUUID()** to request to get the uuid from the device. Register **GetUUIDReceiver** to get the uuid.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetUUID()	Request to get the uuid	GetUUIDReceiver	Get the uuid
(UUID: 0x00002A23-0000-1000-			
8000-00805F9B34FB)			

8. Call the **ISCNIRScanSDK.ReadActivateState()** to request to get whether the device is activate. Register **ReturnReadActivateStatusReceiver** to get whether the device is activated. Call the **ISCNIRScanSDK.SetActiveConfig()** to set the active scan configuration.

(a) Device is activated : Call the **SetDeviceButtonStatus()** to check whether user want to lock button.

(b) Call the ISCNIRScanSDK.ControlPhysicalButton(ISCNIRScanSDK.PhysicalButton.Lock) to lock device button.

(c) Call the ISCNIRScanSDK.ControlPhysicalButton(ISCNIRScanSDK.PhysicalButton.Unlock) to unlock device button.

(d) Device is not activated : Check whether store license in the app. The shared preferences key is **ISCNIRScanSDK.SharedPreferencesKeys.licensekey**.

(1) Have license : Call the ISCNIRScanSDK.SetLicenseKey(data) to set the license. Register RetrunActivateStatusReceiver to get the license is valid or not. Call the SetDeviceButtonStatus() to check whether user want to lock button.

(2) Not have license : Call the **SetDeviceButtonStatus()** to check whether user want to lock button.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.ReadActivateState() (UUID: 0x43484130-444C-5020-4E49- 52204E616E6F)	Request to get whether the device is activate	RetrunReadActivateStatusReceiver	Get whether the device is activated
ISCNIRScanSDK.SetActiveConfig() (UUID: 0x43484118-444C-5020-4E49- 52204E616E6F)	Set the active scan configuration	X	X



SetDeviceButtonStatus()	Check whether user want to lock button	Х	X
ISCNIRScanSDK.ControlPhysicalButton(I SCNIRScanSDK.PhysicalButton.Lock) (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Lock device button	x	X
ISCNIRScanSDK.ControlPhysicalButton(I SCNIRScanSDK.PhysicalButton.Unlock) (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Unlock device button	x	X
ISCNIRScanSDK.SetLicenseKey(data) (UUID: 0x43484130-444C-5020-4E49- 52204E616E6F)	Set the license	RetrunActivateStatusReceiver (UUID: 0x43484131-444C-5020- 4E49-52204E616E6F)	Get the license is valid or not

Perform Scan – Normal

ScanViewActivity.java



- Call the PerformScan(long delaytime) to scan the sample. The delay time is set to 300ms to avoid the BLE hang. In PerformScan(long delaytime), call the ISCNIRScanSDK.StartScan() to trigger the scan event. You should register ScanDataReadyReceiver to get the scan result.
- 2. Call the **ISCNIRScanSDK.GetDeviceStatus()** to request to get the status of the device including battery capacity, system temperature, humidity, total lamp time, the byte of the device status and the byte of the error status. Register **GetDeviceStatusReceiver** to get the device status.
- 3. Call the ISCNIRScanSDK.GetScanLampRampUpADC() to request to get lamp ramp up ADC. Register ReturnLampRampUpADCReceiver to get lamp ramp up ADC.
- 4. Call the **ISCNIRScanSDK.GetLampADCAverage()** to request to get lamp ADC among repeated times. Register **ReturnLampADCAverageReceiver** to get lamp ADC among repeated times.
- 5. Call the writeCSV(Scan_Spectrum_Data) to save the scan report.



API	Process	Receiver	Receiver
			Description
PerformScan(long delaytime)	Set delay time to scan sample to avoid BLE hang.	X	X
ISCNIRScanSDK.StartScan() (UUID: 0x4348411D-444C-5020-4E49- 52204E616E6F 0x43484127-444C-5020-4E49-52204E616E6F)	Trigger the scan event	ScanDataReadyReceiver (UUID: 0x43484128-444C-5020-4E49- 52204E616E6F)	Get the scan result
ISCNIRScanSDK.GetDeviceStatus() (UUID: 0x00002A19-0000-1000-8000-00805F9B34FB 0x43484101-444C-5020-4E49-52204E616E6F 0x43484102-444C-5020-4E49-52204E616E6F 0x43484103-444C-5020-4E49-52204E616E6F 0x43484104-444C-5020-4E49-52204E616E6F)	Request to get the status of the device	GetDeviceStatusReceiver	Get the device status
ISCNIRScanSDK.GetScanLampRampUpADC() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ramp up ADC	ReturnLampRampUpADCReceiver (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Get lamp ramp up ADC
ISCNIRScanSDK.GetLampADCAverage() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ADC among repeated times	ReturnLampADCAverageReceiver (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Get lamp ADC among repeated times
writeCSV(Scan_Spectrum_Data)	Save the scan report	X	Х

Perform Scan – Quick Set

ScanViewActivity.java



1. Set the scan configuration to the device. (Note: This only set the configuration to the device memory and generating the scan patterns accordingly.)



(a) Call the **ChangeScanConfigToByte()** and write quick set UI settings to **ISCNIRScanSDK**. **ScanConfigInfo write_scan_config**. It will return the byte array of the scan config (return **EXTRA_DATA**).

(b) Call the ISCNIRScanSDK.ScanConfig(EXTRA_DATA,ISCNIRScanSDK.ScanConfig.SET) to set the configuration to the device. Register WriteScanConfigStatusReceiver to verify whether the scan configuration was set successfully.

(c) Check whether the set scan configuration in the device memory is valid. Call the **ISCNIRScanSDK.ReadCurrentScanConfig()** to get the current device configuration. Register **ReturnCurrentScanConfigurationDataReceiver** to get the byte array of current configuration of the device.

(d) Call the

Compareconfig(intent.getByteArrayExtra(ISCNIRScanSDK.EXTRA_CURRENT_CONFIG_DATA)) to compare whether the configuration set by the device and the configuration set by the user in the quick set are the same.

- Call the PerformScan(long delaytime) to scan the sample. The delay time is set to 300ms to avoid the BLE hang. In PerformScan(long delaytime), call the ISCNIRScanSDK.StartScan() to trigger the scan event. You should register ScanDataReadyReceiver to get the scan result.
- 3. Call the **ISCNIRScanSDK.GetDeviceStatus()** to request to get the status of the device including battery capacity, system temperature, humidity, total lamp time, the byte of the device status and the byte of the error status. Register **GetDeviceStatusReceiver** to get the device status.
- 4. Call the ISCNIRScanSDK.GetScanLampRampUpADC() to request to get lamp ramp up ADC. Register ReturnLampRampUpADCReceiver to get lamp ramp up ADC.
- 5. Call the ISCNIRScanSDK.GetLampADCAverage() to request to get lamp ADC among repeated times. Register ReturnLampADCAverageReceiver to get lamp ADC among repeated times.
- 6. Call the writeCSV(Scan_Spectrum_Data) to save the scan report.

API	Process	Receiver	Receiver Description
ChangeScanConfigToByte()	Convert scan Config to byte array. Return value is EXTRA_DATA	x	X
ISCNIRScanSDK.ScanConfig(EXTRA_DATA,IS CNIRScanSDK.ScanConfig.SET) (UUID: 0x43484142-444C-5020-4E49- 52204E616E6F)	Set the configuration to the device	WriteScanConfigStatusReceiv er (UUID: 0x43484143-444C- 5020-4E49-52204E616E6F)	Verify whether the scan configuration was set successfully
ISCNIRScanSDK.ReadCurrentScanConfig() (UUID: 0x43484140-444C-5020-4E49- 52204E616E6F)	Get the current device configuration	ReturnCurrentScanConfigurati onDataReceiver (UUID: 0x43484141-444C- 5020-4E49-52204E616E6F)	Get the byte array of current configuration of the device
Compareconfig(intent.getByteArrayExtra(ISCNIR ScanSDK.EXTRA_CURRENT_CONFIG_DATA))	Compare whether the configuration set by the device and the configuration set by the user in the quick set are the same	X	X
PerformScan(long delaytime)	Set delay time to scan sample to avoid BLE hang.	x	X
ISCNIRScanSDK.StartScan() (UUID: 0x4348411D-444C-5020-4E49- 52204E616E6F	Trigger the scan event	ScanDataReadyReceiver (UUID: 0x43484128-444C- 5020-4E49-52204E616E6F)	Get the scan result



0x43484127-444C-5020-4E49-52204E616E6E)			
0,40404121 4440 0020 4240 022042010201)			
ISCNIRScanSDK.GetDeviceStatus()	Request to get the status	GetDeviceStatusReceiver	Get the device status
(UUID:	of the device		
0x00002A19-0000-1000-8000-00805F9B34FB			
0x43484101-444C-5020-4E49-52204E616E6F			
0x43484102-444C-5020-4E49-52204E616E6F			
0x43484109-444C-5020-4E49-52204E616E6F			
0x43484103-444C-5020-4E49-52204E616E6F			
0x43484104-444C-5020-4E49-52204E616E6F)			
ISCNIRScanSDK.GetScanLampRampUpADC()	Request to get lamp ramp	ReturnLampRampUpADCRec	Get lamp ramp up
(UUID: 0x4348410B-444C-5020-4E49-	up ADC	eiver	ADC
52204E616E6F)		(UUID: 0x4348410B-444C-	
,		5020-4E49-52204E616E6F)	
ISCNIRScanSDK.GetLampADCAverage()	Request to get lamp ADC	ReturnLampADCAverageRec	Get lamp ADC among
(UUID: 0x4348410B-444C-5020-4E49-	among repeated times	eiver	repeated times
52204E616E6F)	5,	(UUID: 0x4348410B-444C-	
,		5020-4E49-52204E616E6F)	
writeCSV(Scan_Spectrum_Data)	Save the scan report	Х	Х
-			

Perform Scan – Manual

ScanViewActivity.java



- 1. Turn on the manual mode will support to set lamp on or off, scan PGA and scan repeats.
- 2. Manual Setting



- (a) Call the ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.ON) to turn on the lamp. Register ReturnSetLampReceiver to know the set is completed.
- (b) Call the ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.OFF) to turn off the lamp. Register ReturnSetLampReceiver to know the set is completed.
- (c) Call the **ISCNIRScanSDK.SetPGA(pga)** to set scan PGA. Register **ReturnSetPGAReceiver** to know the set is completed.
- (d) Call the ISCNIRScanSDK.SetScanRepeat(scan_repeat) to set scan repeats. Register ReturnSetScanRepeatsReceiver to know the set is completed.
- 3. Call the **PerformScan(long delaytime)** to scan the sample. The delay time is set to 300ms to avoid the BLE hang. In **PerformScan(long delaytime)**, call the **ISCNIRScanSDK.StartScan()** to trigger the scan event. You should register **ScanDataReadyReceiver** to get the scan result.
- 4. Call the **ISCNIRScanSDK.GetDeviceStatus()** to request to get the status of the device including battery capacity, system temperature, humidity, total lamp time, the byte of the device status and the byte of the error status. Register **GetDeviceStatusReceiver** to get the device status.
- 5. Call the ISCNIRScanSDK.GetScanLampRampUpADC() to request to get lamp ramp up ADC. Register ReturnLampRampUpADCReceiver to get lamp ramp up ADC.
- 6. Call the **ISCNIRScanSDK.GetLampADCAverage()** to request to get lamp ADC among repeated times. Register **ReturnLampADCAverageReceiver** to get lamp ADC among repeated times.
- 7. Call the writeCSV(Scan_Spectrum_Data) to save the scan report.

Close the manual mode as follows or reference to ChangeLampState():

- If the lamp is turned on, should close the lamp(ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.OFF)).
- Change the lamp state to auto(ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.AUTO)).

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK .LampState.ON) (UUID: 0x43484144-444C-5020-4E49- 52204E616E6F)	Turn on the lamp	ReturnSetLampReceiver	Know the set is completed
ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK .LampState.OFF) (UUID: 0x43484144-444C-5020-4E49- 52204E616E6F)	Turn off the lamp	ReturnSetLampReceiver	Know the set is completed
ISCNIRScanSDK.SetPGA(pga) (UUID: 0x43484146-444C-5020-4E49- 52204E616E6F)	Set scan PGA	ReturnSetPGAReceiver	Know the set is completed
ISCNIRScanSDK.SetScanRepeat(scan_repeat) (UUID: 0x43484147-444C-5020-4E49- 52204E616E6F)	Set scan repeats	ReturnSetScanRepeatsReceiv er	Know the set is completed
PerformScan(long delaytime)	Set delay time to scan sample to avoid BLE hang.	X	X
ISCNIRScanSDK.StartScan() (UUID: 0x4348411D-444C-5020-4E49- 52204E616E6F 0x43484127-444C-5020-4E49-52204E616E6F)	Trigger the scan event	ScanDataReadyReceiver (UUID: 0x43484128-444C- 5020-4E49-52204E616E6F)	Get the scan result
ISCNIRScanSDK.GetDeviceStatus() (UUID: 0x00002A19-0000-1000-8000-00805F9B34FB 0x43484101-444C-5020-4E49-52204E616E6F 0x43484102-444C-5020-4E49-52204E616E6F	Request to get the status of the device	GetDeviceStatusReceiver	Get the device status



0x43484109-444C-5020-4E49-52204E616E6F			
0x43484103-444C-5020-4E49-52204E616E6F			
0x43484104-444C-5020-4E49-52204E616E6F)			
ISCNIRScanSDK.GetScanLampRampUpADC() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ramp up ADC	ReturnLampRampUpADCRec eiver (UUID: 0x4348410B-444C- 5020-4E49-52204E616E6F)	Get lamp ramp up ADC
ISCNIRScanSDK.GetLampADCAverage() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ADC among repeated times	ReturnLampADCAverageRec eiver (UUID: 0x4348410B-444C- 5020-4E49-52204E616E6F)	Get lamp ADC among repeated times
writeCSV(Scan_Spectrum_Data)	Save the scan report	X	Х
ChangeLampState()	Leaving the manual scan page to restore the lamp state to auto	x	x

Update Reference to the Device

ScanViewActivity.java



- 1. Set reference scan configuration.
 - (a) Call the ISCNIRScanSDK.SetReferenceParameter(int MINWAV,int MAXWAV) to set the reference configuration to the device. Register WriteScanConfigStatusReceiver to get back and check whether the set scan configuration in the device memory is valid.

Note: In our APP, Need to set the reference set flag to Compareconfig(byte EXTRA_DATA[]) to determine that the reference config needs to be compared now.

reference_set_config = true;



(b) Call the ISCNIRScanSDK.ReadCurrentScanConfig() to get the current device configuration. Register ReturnCurrentScanConfigurationDataReceiver to get the byte array of current configuration of the device.

(c) Call the

Compareconfig(intent.getByteArrayExtra(ISCNIRScanSDK.EXTRA_CURRENT_CONFIG_DAT A)) to compare whether the configuration set by the device and the configuration set by the reference are the same.

- Call the PerformScan(long delaytime) to scan the reference sample. The delay time is set to 300ms to avoid the BLE hang. In PerformScan(long delaytime), call the ISCNIRScanSDK.StartScan() to trigger the scan event. You should register ScanDataReadyReceiver to get the scan results.
- 3. Call the **ISCNIRScanSDK.GetDeviceStatus()** to request to get the status of the device including battery capacity, system temperature, humidity, total lamp time, the byte of the device status and the byte of the error status. Register **GetDeviceStatusReceiver** to get the device status.
- 4. Call the ISCNIRScanSDK.GetScanLampRampUpADC() to request to get lamp ramp up ADC. Register ReturnLampRampUpADCReceiver to get lamp ramp up ADC.
- 5. Call the **ISCNIRScanSDK.GetLampADCAverage()** to request to get lamp ADC among repeated times. Register **ReturnLampADCAverageReceiver** to get lamp ADC among repeated times.
- 6. Call the **writeCSV(Scan_Spectrum_Data)** to save the scan report. In our APP, will change configuration name to Reference.
- 7. Call the ISCNIRScanSDK.SaveReference() to set the reference to the device.

Note: In our APP, will pop out the confirm dialog and should set reference flag.

saveReference = true;

 After finish saving the reference to the device, call the ISCNIRScanSDK.ScanConfig(ActiveConfigByte,ISCNIRScanSDK.ScanConfig.SET) to set active scan configuration to the device and disconnect.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.SetReferenceParameter(int MINWAV,int MAXWAV) (UUID: 0x43484142-444C-5020-4E49- 52004E645E6E	Set the reference configuration to	WriteScanConfigStatusReceiver (UUID: 0x43484142-444C-5020- 4E49-52204E616E6F)	Get back and check whether the set scan configuration in the
S2204E616E6F) ISCNIRScanSDK.ReadCurrentScanConfig() (UUID: 0x43484140-444C-5020-4E49- 52204E616E6F)	Get the current device configuration	ReturnCurrentScanConfigurationData Receiver (UUID: 0x43484141-444C-5020- 4E49-52204E616E6F)	Get the byte array of current configuration of the device
Compareconfig(intent.getByteArrayExtra(ISCNIRSc anSDK.EXTRA_CURRENT_CONFIG_DATA))	Compare whether the configuration set by the device and the configuration set by the reference are the same	x	x
PerformScan(long delaytime)	Set delay time to scan sample to avoid BLE hang.	X	X



ISCNIRScanSDK.StartScan() (UUID: 0x4348411D-444C-5020-4E49-	Trigger the scan event	ScanDataReadyReceiver (UUID: 0x43484128-444C-5020-	Get the scan result
0x43484127-444C-5020-4E49-52204E616E6F)		4E49-52204E616E6F)	
ISCNIRScanSDK.GetDeviceStatus() (UUID: 0x00002A19-0000-1000-8000-00805F9B34FB 0x43484101-444C-5020-4E49-52204E616E6F 0x43484102-444C-5020-4E49-52204E616E6F 0x43484103-444C-5020-4E49-52204E616E6F 0x43484104-444C-5020-4E49-52204E616E6F	Request to get the status of the device	GetDeviceStatusReceiver	Get the device status
ISCNIRScanSDK.GetScanLampRampUpADC() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ramp up ADC	ReturnLampRampUpADCReceiver (UUID: 0x4348410B-444C-5020- 4E49-52204E616E6F)	Get lamp ramp up ADC
ISCNIRScanSDK.GetLampADCAverage() (UUID: 0x4348410B-444C-5020-4E49- 52204E616E6F)	Request to get lamp ADC among repeated times	ReturnLampADCAverageReceiver (UUID: 0x4348410B-444C-5020- 4E49-52204E616E6F)	Get lamp ADC among repeated times
writeCSV(Scan_Spectrum_Data)	Save the scan report	X	x
ISCNIRScanSDK.SaveReference() (UUID: 0x43484132-444C-5020-4E49- 52204E616E6F)	Set the reference to the device	x	X
ISCNIRScanSDK.ScanConfig(ActiveConfigByte,IS CNIRScanSDK.ScanConfig.SET) (UUID: 0x43484142-444C-5020-4E49- 52204E616E6F)	Set active scan configuration to the device	x	X

Activation Key Setting

ActivationViewActivity.java



1. Call the ISCNIRScanSDK.SetLicenseKey(data) to set the activation key(key length is 24). Register ReturnActivateStatusReceiver to check whether the activation key is set success.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.SetLicenseKey(d ata) (UUID: 0x43484130-444C-5020- 4E49-52204E616E6F)	Set the activation key	ReturnActivateStatusReceiver (UUID: 0x43484131-444C-5020- 4E49-52204E616E6F)	Check whether the activation key is set success



Device Information

DeviceInfoViewActivity.java



 Call the ISCNIRScanSDK.GetDeviceInfo() to request to get the device information including the device of manufacturer, model name, serial number, HW version and TIVA version. Register DeviceInfoReceiver to get the device information.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetDeviceInfo()	Request to get the	DeviceInfoReceiver	Get the device information
(UUID: 0x00002A29-0000-1000-8000-	device information		
00805F9B34FB			
0x00002A24-0000-1000-8000-00805F9B34FB			
0x00002A25-0000-1000-8000-00805F9B34FB			
0x00002A27-0000-1000-8000-00805F9B34FB			
0x00002A26-0000-1000-8000-00805F9B34FB		7	
0x00002A28-0000-1000-8000-00805F9B34FB)			

Device Status

DeviceStatusViewActivity.java



1. Call the **ISCNIRScanSDK.GetDeviceStatus()** to request to get the status of the device including battery capacity, system temperature, humidity, total lamp time, the byte of the device status and the byte of the error status. Register **mStatusReceiver** to get the device status.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetDeviceStatus()	Request to get the	mStatusReceiver	Get the device status
(UUID:	status of the device		
0x00002A19-0000-1000-8000-00805F9B34FB			
0x43484101-444C-5020-4E49-52204E616E6F			
0x43484102-444C-5020-4E49-52204E616E6F			
0x43484109-444C-5020-4E49-52204E616E6F			
0x43484103-444C-5020-4E49-52204E616E6F			
0x43484104-444C-5020-4E49-52204E616E6F)			



List Scan Configuration

ScanConfigurationsViewActivity.java



- Call the ISCNIRScanSDK.GetScanConfig() to request to get the number of scan configuration and scan configuration data. Register ScanConfSizeReceiver to get the number of scan configuration. Register ScanConfReceiver to get scan configuration data.
- 2. Call the **ISCNIRScanSDK.GetActiveConfig()** to request to get the active configuration index in the device. Register **GetActiveScanConfReceiver** to get active configuration index.

API	Process	Receiver	Receiver Description
ISCNIRScanSDK.GetScanConfig() (UUID:	Get the number of scan configuration->Get scan	ScanConfSizeReceiver	Get the number of scan configuration
0x43484113-444C-5020-4E49- 52204E616E6F 0x43484114-444C-5020-4E49- 52204E616E6F)	configuration data	ScanConfReceiver (UUID: 0x43484115-444C-5020- 4E49-52204E616E6F)	Get scan configuration data
ISCNIRScanSDK.GetActiveConfig() (UUID: 0x43484118-444C-5020- 4E49-52204E616E6F)	Request to get the active configuration index	GetActiveScanConfReceiver (UUID: 0x43484118-444C- 5020-4E49-52204E616E6F)	Get active configuration index

Add Scan Configuration to the Device

AddScanConfigViewActivity.java



 Call the ChangeScanConfigToByte() and write UI settings to ISCNIRScanSDK.WriteScanConfiguration(write_scan_config). It will return the byte array of the scan configuration (return EXTRA_DATA).



 Call the ISCNIRScanSDK.ScanConfig(EXTRA_DATA,ISCNIRScanSDK.ScanConfig.SAVE) to set the configuration to the device. Register WriteScanConfigStatusReceiver to check whether the setting is successful.

API	Process	Receiver	Receiver Description
ChangeScanConfigToByte()	Transmit scan configuration to the byte that user set.	X	X
ISCNIRScanSDK.ScanConfig(EXT RA_DATA,ISCNIRScanSDK.Scan Config.SAVE) (UUID: 0x43484142-444C-5020- 4E49-52204E616E6F)	Set the configuration to the device	WriteScanConfigStatusReceiver (UUID: 0x43484143-444C-5020- 4E49-52204E616E6F)	Check whether the setting is successful

Warm-Up the Device

HomeViewActivity.java

When the user clicks on the connection, the Intent will be used to pass the flag whether the user wants to warm up the device. See **newscanhIntent.putExtra("warmup",switch_Warmup.isChecked())**;. After clicks on the connection, see the section of <u>Connected to the Device 1-(b)</u>. It will teach how to open the lamp to do the warm up the device.

ScanViewActivity.java

Warm-up time reaches three minutes, you can directly switch pages or scan, and the device will automatically turn off the lamp. See **ChangeLampState()**. If device have do the warm up, will call **ISCNIRScanSDK.ControlLamp(ISCNIRScanSDK.LampState.AUTO)** to close the lamp.

Warning: When the user sets connect with lamp on, must pay attention to the warm-up time and confirm that the device is warmed up and turn off the lamp before leaving.



Calculate AU from a New Reference Scan



- Place the reference sample and call the PerformScan (long delaytime) to scan the reference sample. The delay time is set to 300ms to avoid the BLE hang. In PerformScan(long delaytime), call the ISCNIRScanSDK.StartScan() to trigger the scan event. You should register ScanDataReadyReceiver to get the scan result.
- 2. The acquired spectral data(intensity) needs to be temporarily stored in ReferenceIntensity.

For example :

}

```
public static List<Integer> ReferenceIntensity = new ArrayList<Integer>();
```

ReferenceIntensity.clear(); Scan_Spectrum_Data = new ISCNIRScanSDK.ScanResults(Interpret_wavelength,Interpret_intensity,Interpret_uncalibratedIntens ity,Interpret_length);

for (index = 0; index < Scan_Spectrum_Data.getLength(); index++) {
 ReferenceIntensity.add(Scan_Spectrum_Data.getUncalibratedIntensity()[index]);</pre>

 Place the sample and call the PerformScan(long delaytime) to scan the sample. The delay time is set to 300ms to avoid the BLE hang. In PerformScan(long delaytime), call the ISCNIRScanSDK.StartScan() to trigger the scan event. You should register ScanDataReadyReceiver to get the scan result.



4. Use the stored **new reference sample signal intensity and sample signal intensity** to calculate absorbance and reflectance of the sample.

For example :

Scan_Spectrum_Data = new

ISCNIRScanSDK.ScanResults(Interpret_wavelength,Interpret_intensity,Interpret_uncalibratedIntensit y,Interpret_length);

mXValues.clear(); mIntensityFloat.clear(); mAbsorbanceFloat.clear(); mReflectanceFloat.clear(); mWavelengthFloat.clear(); mReferenceFloat.clear(); int index; for (index = 0; index < Scan_Spectrum_Data.getLength(); index++) { mXValues.add(String.format("%.02f", ISCNIRScanSDK.ScanResults.getSpatialFreq(mContext, Scan_Spectrum_Data.getWavelength()[index])));

mIntensityFloat.add(new Entry((float) Scan_Spectrum_Data.getWavelength()[index],(float) Scan_Spectrum_Data.getUncalibratedIntensity()[index]));

mAbsorbanceFloat.add(new Entry((float) Scan_Spectrum_Data.getWavelength()[index],(-1) * (float) Math.log10((double) Scan_Spectrum_Data.getUncalibratedIntensity()[index] / (double) ReferenceIntensity.get(index))));

mReflectanceFloat.add(new Entry((float) Scan_Spectrum_Data.getWavelength()[index],(float) Scan_Spectrum_Data.getUncalibratedIntensity()[index] / (float) ReferenceIntensity.get(index)));

mWavelengthFloat.add((float) Scan_Spectrum_Data.getWavelength()[index]);

mReferenceFloat.add(new Entry((float) Scan_Spectrum_Data.getWavelength()[index],(float) **ReferenceIntensity.get(index)**));

Note: The ScanDataReadyReceiver needs to judge whether the scan is to scan the reference sample or the sample to decide whether to go step2 or step4. In addition, the ReferenceIntensity can also be saved to the local file, and the local file can be read and loaded into the ReferenceIntensity before entering the sample scanning. In this demonstration, it should be noted that the scan configuration (including PGA) of the scanned reference sample must be the same as the scan configuration of the current scanned sample.